

Anti-human TAU MTBR mab 1F8

Catalogue number	Package size
847-0102010901	100 µg
847-0102010903	1 mg

For research use only.

Product:	monoclonal antibody
Clone:	1F8
Host:	mouse
Isotype:	IgG1 _k
Immunogen:	human TAU sequence K12CF3R
Specificity:	SKVTSKCG of human TAU (2N3R) amino acids 317-323
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	Carbonate buffer pH 9.6 without additives
Storage:	at -20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot, Immunohistochemistry
Dilution:	1 µg/ml in Western Blot
Reactivity:	Microtubuli binding region TAU protein



Roboscreen GmbH
Hohmannstraße 7
04129 Leipzig • Germany

Phone: +49 341 989 734 0
Fax: +49 341 989 734 199
info@roboscreen.com
www.roboscreen.com

CEO
Dr. Ingolf Lachmann

Ust-IdNr. DE 813064314
County Court Leipzig
HRB 17503
Leipzig